

Technical Specification

PRODUCT : SODIUM HYPOCHLORITE SOLUTION IN WATER RE - Pure**CODE : 370320****MOLECULAR WEIGHT :74,44****CAS N° : 7681-52-9****FORMULA : NaClO**

TEST	U.M.	SPECIFICATION
Description	-	Yellow clear liquid
Alcalinity (NaOH)	% (w/w)	<= 1.8
Assay (iodometric)	% (Cl)	5 - 9
Density at 20°C	-	1.1 - 1.2

The content stated refers to the filling day. Depending on storage period and temperature a decrease in content is to be expected.
Store at ambient temperature. Keep away from heat.

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in