

Technical Specification

PRODUCT : TOLUENE RE - Pure**CODE : 386000****MOLECULAR WEIGHT :92,14****CAS N° : 108-88-3****FORMULA : C₆H₅CH₃**

TEST	U.M.	SPECIFICATION
Description	-	Clear, colourless liquid
Identification	-	Positive
Density d ₂₀ /4	-	0.862 - 0.872
Refractive index at 20°C	-	1.491 - 1.501
Boiling point	°C	109.9 - 111.4
Water (K.F.)	ppm	<= 300
Residue on evaporation	ppm	<= 50
Benzene	%	<= 0.02
Total sulphur (S)	ppm	<= 100
Assay (GC)	%	>= 99.5

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in