

# Technical Specification

**PRODUCT : TOLUENE RE - ASTM****CODE : 386100****MOLECULAR WEIGHT :92.14****CAS N° : 108-88-3****FORMULA : C<sub>6</sub>H<sub>5</sub>CH<sub>3</sub>**

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless to pale yellow liquid
Peroxide	ppm	<= 5
Density at 15°C	-	-
Water (K.F.)	ppm	<= 200
Toluene (ASTM)	% (V/V)	>= 99.5
Antioxidant	mg/L	1500 - 2500

Suitable for ASTM methods D2700 and D2699.

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.  
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in