

Technical Specification

PRODUCT : Xylene, mix of isomers RS For histology**CODE : 392620****MOLECULAR WEIGHT :106,17****CAS N° : 1330-20-7****FORMULA : C8H10**

TEST	U.M.	SPECIFICATION
------	------	---------------

Description	-	Clear liquid
Identification (I.R.)	-	Positive
Benzene	ppm	<= 50
Assay (isomeric mixture)	%	>= 98.5

General description :
Boiling point : 140 °C

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in