

Technical Specification

PRODUCT : SULFURIC ACID 90% RS - For analysis according to Gerber**CODE : 410390****MOLECULAR WEIGHT :98,08****CAS N° : 7664-93-9****FORMULA : H2SO4**

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Density at 20°C	-	1.815 - 1.825
Assay (acidimetric)	%	90 - 92

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in