

Technical Specification

PRODUCT : TRIFLUOROACETIC ACID RS - For LC/MS**CODE : 411540**

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Identification	-	Positive
Water (K.F.)	%	<= 0.05
Chloride	ppm	<= 10
Fluoride	ppm	<= 50
Sulfate	ppm	<= 10
Assay (acidimetric)	%	>= 99.9
Suitability	-	LC-MS Tested

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in