

Technical Specification

PRODUCT : ACETONITRILE RS - For preparative HPLC

CODE : 412400

MOLECULAR WEIGHT :41.05

CAS N° : 75-05-8

FORMULA : CH₃CN

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Identification	-	Positive
Colour	APHA	<= 10
Density at 20°C	-	0.781 - 0.785
Refractive index at 20°C	-	1.3420 - 1.3440
Boiling point	°C	81.1 - 82.1
Acidity or alkalinity	meq/g	<= 0.0008
Water (K.F.)	ppm	<= 300
Residue on evaporation	ppm	<= 5
Assay (GC)	%	>= 99.9
UV Transmittance	-	-
at 250 nm	%	>= 98
at 230 nm	%	>= 50

Filtered at 1 micron

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
 CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in