

Technical Specification

PRODUCT : Propan-2-ol RS MOS For electronic use

CODE : 415160

MOLECULAR WEIGHT :60,1

CAS N° : 67-63-0

FORMULA : CH₃CHOHCH₃

TEST	U.M.	SPECIFICATION
Description	-	Clear liquid
Colour	APHA	<=10
Identification (I.R.)	-	Positive
Water miscibility	-	Conform
Resistivity	MOhm.cm	>= 10
Assay (GC)	%	>= 99.7
Density at 20°C	-	0.784 - 0.786
Boiling point	°C	82.1 - 82.6
Water (K.F.)	ppm	<= 500
Residue on evaporation	ppm	<= 5
Acidity (propionic acid)	ppm	<= 10
Alcalinity (NH ₃)	ppm	<= 2
Aldehydes - ketones	ppm	<= 50
Chloride	ppm	<= 0.2
Phosphate	ppm	<= 0.5
Heavy metals (as Pb)	ppm	<= 0.1
Subs. reducing KMnO ₄	ppm	<= 2.5
Total sulphur	ppm	<= 1
Ag	ppm	<= 0.02
Al	ppm	<= 0.05
As	ppm	<= 0.01
Au	ppm	<= 0.05
B	ppm	<= 0.01
Ba	ppm	<= 0.1
Be	ppm	<= 0.02
Bi	ppm	<= 0.02
Ca	ppm	<= 0.2
Cd	ppm	<= 0.01
Co	ppm	<= 0.01
Cr	ppm	<= 0.01
Cu	ppm	<= 0.01
Fe	ppm	<= 0.1
Ga	ppm	<= 0.02
In	ppm	<= 0.02
K	ppm	<= 0.1
Li	ppm	<= 0.02
Mg	ppm	<= 0.1
Mn	ppm	<= 0.01
Mo	ppm	<= 0.05
Na	ppm	<= 0.5
Ni	ppm	<= 0.01
Pb	ppm	<= 0.01
Pt	ppm	<= 0.05
Sb	ppm	<= 0.01
Sn	ppm	<= 0.02

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.

CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in

Technical Specification

TEST	U.M.	SPECIFICATION
Sr	ppm	<= 0.02
Ti	ppm	<= 0.05
Tl	ppm	<= 0.05
V	ppm	<= 0.05
Zn	ppm	<= 0.01
Zr	ppm	<= 0.05
Filtered to 0.2 micron		

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
 CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in