

# Technical Specification

**PRODUCT : SILVER NITRATE 0.5 N (0.5 mol/L) RPE - For analysis****CODE : 424050000**

| TEST                  | U.M. | SPECIFICATION           |
|-----------------------|------|-------------------------|
| Description           | -    | Clear colourless liquid |
| Assay (potentiometry) | N    | 0.4995 - 0.5005         |

This volumetric solution is checked by potassium chloride volumetric standard, traceable to the SRM 999c of NIST (National Institute of Standards and Technology)

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.  
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in