

Technical Specification

**PRODUCT : L(+)
Arginine RPE - For analysis****CODE : 424270****MOLECULAR WEIGHT :174,20****CAS N° : 74-79-3****FORMULA : C6H14N4O2**

TEST	U.M.	SPECIFICATION
Description	-	White crystalline powder
Identification	-	Positive
Specific optical rotation	°	+26.3 to +27.7
Assay (on dry)	%	>= 98

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in