

# Technical Specification

**PRODUCT : Ehrlich's reagent RS - For microscopy**

**CODE : 446300**

TEST	U.M.	SPECIFICATION
Description	-	clear slightly yellow liquid
Max absorption lambda max (H2O)	nm	454 - 458
ABS (Lambda max; 2:100 H2O)	-	0.37 - 0.41

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.  
 CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in