

Technical Specification

PRODUCT : May Grünwald reagent RS - For hematology

CODE : 460580

TEST	U.M.	SPECIFICATION
Description	-	Blue clear liquid
Max lambda abs. max1 (MetOH;dil.1,25:250)	nm	648 - 653
Lambda abs. max1 (MetOH;dil.1,25:250)	-	1.517 - 1.862
Max lambda abs. Max2 (MetOH;dil.1,25:250)	nm	520 - 530
Lambda abs. max2 (MetOH;dil.1,25:250)	-	0.827 - 1.103

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
 CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in