

Technical Specification

PRODUCT : DICHLOROMETHANE RS for GC-MS ,stabilized with ethanol

CODE : 463330

MOLECULAR WEIGHT :84.93

CAS N° : 75-09-2

FORMULA : CH₂Cl₂

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Refractive index at 20°C	-	1.422 - 1.426
Water (K.F.)	ppm	<= 100
Residue on evaporation	ppm	<= 3
Acidity (HCl)	ppm	<= 5
Colour	APHA	<= 10
Assay (GC) (w ithout stabilizer)	%	>= 99.95
Stabilizer (Ethanol)	% (w /w)	0.1 - 0.4
Carbon tetrachloride	ppm	<= 20
Chloroform	ppm	<= 40
GC-MS.Individ. peak (hexadecane)	µg/L	<= 2
Ret. range n-undecane to n-tetracontane	-	-
(scanning area 30-600 amu)	-	-

Suitable for residue analysis.

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
 CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in