

Technical Specification

PRODUCT : PYRIDINE RS - For titration according to Karl Fischer**CODE : 469650****MOLECULAR WEIGHT :79.10****CAS N° : 110-86-1****FORMULA : C₅H₅N**

| TEST | U.M. | SPECIFICATION |
|----------------------------|------|-------------------------|
| Description | - | Clear colourless liquid |
| Identification | - | Positive |
| Density d ₂₀ /4 | - | 0.979 - 0.985 |
| Refractive index at 20°C | - | 1.5050 - 1.5140 |
| Boiling point | °C | 114.2 - 116.2 |
| Water (K.F.) | ppm | <= 500 |
| Assay (GC) | % | >= 99.6 |

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in