

Technical Specification

PRODUCT : Potassium acetate RPE - For analysis

CODE : 470140

MOLECULAR WEIGHT :98,15

CAS N° : 127-08-2

FORMULA : CH₃COOK

TEST	U.M.	SPECIFICATION
Description	-	White granular pow der
Identification (I.R.)	-	Positive
pH sol. 5% at 20°C	-	7.5 - 8.5
Chloride	ppm	<= 50
Sulfate	ppm	<= 50
Heavy metals (as Pb)	ppm	<= 10
Fe	ppm	<= 20
Zn	ppm	<= 20
Assay (on dry)	%	>= 99
Water content	%	<= 1.0
Ca	ppm	<= 100
Na	%	<= 0.4
Mg	ppm	<= 100
Hg	ppm	<= 1

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in