

Technical Specification

PRODUCT : POTASSIUM DICHROMATE 0.25 N (0.0417 mol/L) RS - For environmental analysis (COD determination)

CODE : 470450

TEST	U.M.	SPECIFICATION
Description	-	Orange clear liquid
Assay (potentiometry)	N	0.2495 - 0.2505

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in