

# Technical Specification

**PRODUCT : POTASSIUM PERMANGANATE 1N (0.2 mol/L) RPE - For analysis**

**CODE : 473510000**

TEST	U.M.	SPECIFICATION
Description	-	Clear purple liquid
Assay (potentiometry)	N	0.99 - 1.01

This volumetric solution is checked by sodium thiosulfate solution, tested with potassium dichromate volumetric standard traceable to the SRM 136f of NIST (National Institute of Standards and Technology)

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.  
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in