

Technical Specification

PRODUCT : Safranine T hydroalcoholic solution for Gram-Hucker Kit RS - For bacteriology

CODE : 477240

TEST	U.M.	SPECIFICATION
Description	-	Red clear liquid
Maximum lambda absorption (ETOH 50%)	nm	530 - 534
Abs.(Lambda max;ETOH50%;0,003g/L)	-	0.378 - 0.471

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
 CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in