

# Technical Specification

**PRODUCT : Redox solution 468 mV at 25°C RS For electrochemistry**

**CODE : 478050**

TEST	U.M.	SPECIFICATION
Description	-	Clear green light brown liquid
Identification	-	Positive
Redox potential at 25°C	mV	458 - 478

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.  
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in