

# Technical Specification

**PRODUCT : BUFFER SOLUTION pH 3 RPE - For analysis**

**CODE : 486250**

TEST	U.M.	SPECIFICATION
Clear, colourless solution	-	Conform
pH	unite pH	2.98 - 3.02
Temperature of measurement	°C	19 - 21

Accuracy : +/- 0.02 pH units - Standardized against NIST reference materials.

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.  
 CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in