

Technical Specification

PRODUCT : BUFFER SOLUTION pH 6 RPE - For analysis**CODE : 486330**

| TEST | U.M. | SPECIFICATION |
|----------------------------|----------|---------------|
| Clear, colourless solution | - | Conform |
| pH | unite pH | 5.98 - 6.02 |
| Temperature of measurement | °C | 19 - 21 |

Accuracy : +/- 0.02 pH units - Standardized against NIST reference materials.

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in