

# Technical Specification

**PRODUCT : Mixture for checking solderings RS - For agroalimentary analysis****CODE : 502670**

TEST	U.M.	SPECIFICATION
Description	-	Clear purple liquid
Density at 20°C	-	0.788 - 0.798

Sodium sulfosuccinate 65% : 32 g

Fuschine basique : 4 g

Isopropanol : qsp 1000 mL

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.  
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in