

Technical Specification

PRODUCT : HEPTANE mixture of isomers RE - Pure**CODE : 524411**

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Colour	APHA	<= 10
Density d20/4	-	1.070 - 1.080
Assay (GC)	-	Conform
Composition:		
Chloroform stab Amylène	: 490 mL	
n Heptane	: 490 mL	
Methyl alcohol	: 20 mL	

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in