

# Technical Specification

**PRODUCT : OCTANE 80 Blend RE ASTM****CODE : 525990****MOLECULAR WEIGHT :114.23****CAS N° : 111-65-9****FORMULA : CH<sub>3</sub>(CH<sub>2</sub>)<sub>6</sub>CH<sub>3</sub>**

TEST	U.M.	SPECIFICATION
Description	-	Clear, colourless liquid
n-Heptane (ASTM)	% (V/V)	19.9 - 20.1
Isooctane (ASTM)	% (V/V)	79.9 - 80.1

Suitable for ASTM methods D2700 and D2699.

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.  
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in