

Technical Specification

PRODUCT : Xylene RE Pure Low content in benzene

CODE : 528250

MOLECULAR WEIGHT :106,17

CAS N° : 1330-20-7

FORMULA : C₈H₁₀

TEST	U.M.	SPECIFICATION
Description	-	Clear liquid
Identification (I.R.)	-	Positive
Density at 20°C	-	0.862 - 0.872
Refractive index at 20°C	-	1.4917 - 1.5017
Residue on evaporation	ppm	<= 100
Water (K.F.)	ppm	<= 200
Benzene	ppm	<= 50
Total sulphur	ppm	<= 100
Assay (isomeric mixture)	%	>= 98.5

General description :

Boiling point : 140 °C

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
 CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in