

# Technical Specification

**PRODUCT : Dichloromethane RE Pure stabilized with amylene**

**CODE : 528460**

**MOLECULAR WEIGHT :84,93**

**CAS N° : 75-09-2**

**FORMULA : CH<sub>2</sub>Cl<sub>2</sub>**

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Colour	APHA	<= 10
Density d <sub>20</sub> /4	-	1.320 - 1.330
Refractive index at 20°C	-	1.4194 - 1.4294
Boiling point	°C	39.1 - 40.6
Water (K.F.)	ppm	<= 200
Acidity	ppm	<= 10
Residue on evaporation	ppm	<= 20
Amylene	ppm	20 - 60
Assay (GC)	%	>= 99.8

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.  
 CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in