

Technical Specification

PRODUCT : Standard solution 10000 µS/cm RS - For conductivity

CODE : 575120

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Conductivity at 25°C	µS/cm	9900 - 10100

Traceable to Standard NIST N° 999a
 Confidence interval: ± 20 µS /cm

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
 CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in