

# Technical Specification

**PRODUCT : Standard solution 450000 µS/cm RS - For conductivity**

**CODE : 575210**

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Identification	-	
Conductivity at 25°C	µS/cm	495000 - 505000

Traceable to Standard NIST N° 999a  
 Confidence interval: ± 1000 µS/cm

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.  
 CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in