

Technical Specification

PRODUCT : ACETIC ACID 30% RS - For analysis according to Ph. Eur. Chap. 4.1.1

CODE : 611000401

MOLECULAR WEIGHT :60.05

CAS N° : 64-19-7

FORMULA : CH₃COOH

TEST	U.M.	SPECIFICATION
------	------	---------------

Assay	g/L	290 - 310
-------	-----	-----------

The reagent is specially intended for use in testing Ph Eur and monographs. This solution was produced according 4.1.1 of Ph Eur specifications using Ph Eur reagent quality components and water meeting the Ph Eur monograph requirement for Purified Water.

Instruction for use

The solution bottle should be open for the minimum time required to dispense the solution.

After use, the bottle should be tightly recapped and stored under normal laboratory conditions.

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in