

Compatibility of LLG-Pipette tips "economy" with Eppendorf Pipettes

Pipettor		Eppendorf																																							
Brand		Reference														Research				Research plus																					
Model		variable/single channel						fix/single channel								variable/single channel				variable/single channel																					
Volume Range		0.1-2.5 µl (anthracit)	0.5-10 µl (gray)	2-20 µl (gray)	2-20 µl (yellow)	10-100 µl (yellow)	50-200 µl (yellow)	100-1000 µl (blue)	500-2500 µl (red)	1 µl (gray)	2 µl (gray)	5 µl (gray)	10 µl (gray)	10 µl (yellow)	20 µl (yellow)	25 µl (yellow)	50 µl (yellow)	100 µl (yellow)	200 µl (blue)	250 µl (blue)	500 µl (blue)	1000 µl (blue)	1500 µl (red)	2000 µl (red)	2500 µl (red)	0.1-2.5 µl (dark gray)	0.5-10 µl (gray)	2-20 µl (yellow)	10-100 µl (yellow)	20-200 µl (yellow)	100-1000 µl (blue)	500-5000 µl (violet)	2-20 µl (yellow)	10-100 µl (yellow)	20-200 µl (yellow)	100-1000 µl (blue)					
		Pipette tips "economy" line																																							
Volume	Cat-No.																																								
0.1-10 µl	9409099	■	■	■					■	■	■	■														■										■					
	9409045	■	■	■					■	■	■	■														■										■					
1-200 µl	9409046				■	■	■						■	■	■	■	■	■										■	■	■							■	■			
	6253577				■	■	■						■	■	■	■	■	■										■	■	■								■	■		
	9409047				■	■	■						■	■	■	■	■	■										■	■	■								■	■		
100-1000 µl long	9409048						■												■	■	■	■													■				■		
	9409049						■												■	■	■	■													■				■		
1-5 ml	9409097							■																■		■		■							■						

Compatibility of LLG-Pipette tips "economy" with Gilson Pipettes

Pipettor		Gilson																				
Brand		Pipetman Classic					Fixed Pipetman															
		fix/single channel																				
Volume Range		0.2-2 µl (P2)	1-10 µl (P10)	2-20 µl (P20)	20-100 µl (P100)	50-200 µl (P200)	200-1000 µl (P1000)	1000-5000 µl (P5000)	1-10 ml (P10 ml)	2 µl (F2)	5 µl (F5)	10 µl (F10)	20 µl (F20)	25 µl (F25)	50 µl (F50)	100 µl (F100)	200 µl (F200)	250 µl (F250)	300 µl (F300)	400 µl (F400)	500 µl (F500)	1000 µl (F1000)
Pipette tips "economy" line																						
Volume	Cat-No.																					
0.1-10 µl	9409099		■							■	■	■										
	9409045		■							■	■	■										
1-200 µl	9409046			■	■	■						■		■	■	■						
	6253577			■	■	■						■		■	■	■						
	9409047			■	■	■						■		■	■	■						
100-1000 µl long	9409048						■										■	■	■	■	■	■
	9409049						■										■	■	■	■	■	■
1-10 ml	9409098							■														

Compatibility of LLG-Pipette tips "economy" with Biohit Pipettes

Pipettor		Biohit																															
Brand		Proline																															
Model		variable/single channel						fix/single channel						variable/ multi channel (8)		variable/ multi channel (12)		variable/ electronic single channel		*1	*2												
Volume Range		0.5-10 µl	2-20 µl	5-50 µl	10-100 µl	20-200 µl	100-1000 µl	1000-5000 µl	5 µl	10 µl	20 µl	25 µl	50 µl	100 µl	200 µl	250 µl	500 µl	1000 µl	2000 µl	5000 µl	0.5-10 µl	5-50 µl	50-300 µl	0.5-10 µl	5-50 µl	50-300 µl	5-100 µl	10-500 µl	50-1000 µl	100-5000 µl	5-100 µl	5-100 µl	
Pipette tips "economy" line																																	
Volume	Cat-No.																																
0.1-10 µl	9409099	■							■	■											■			■									
	9409045	■							■	■												■			■								
1-200 µl	9409046		■	■	■	■					■	■	■	■	■							■				■					■	■	
	6253577		■	■	■	■					■	■	■	■	■								■				■				■	■	
	9409047		■	■	■	■					■	■	■	■	■								■				■				■	■	
100-1000 µl long	9409048						■								■	■	■											■					
	9409049						■								■	■	■										■						
1-5 ml	9409097						■												■	■										■			

Compatibility of LLG-Pipette tips "economy" with Brand Pipettes

Pipettor		Brand																										
Brand		Transferpette																										
Model		variable/single channel				fix/single channel						variable/multi channel (8)				variable/multi channel (12)												
Volume Range		0.1-1 µl	0.5-10 µl	2-20 µl	5-50 µl	10-100 µl	100-1000 µl	500-5000 µl	10 µl	50 µl	100 µl	200 µl	200 µl	250 µl	500 µl	1000 µl	2000 µl	0.5-10 µl	5-50 µl	10-100 µl	20-200 µl	30-300 µl	2-20 µl	10-100 µl	20-200 µl	30-300 µl		
Pipette tips "economy" line																												
Volume	Cat-No.																											
0.1-10 µl	9409099	■	■						■									■										
	9409045	■	■						■									■										
1-200 µl	9409046			■	■	■				■	■	■								■	■	■				■	■	
	6253577			■	■	■				■	■	■								■	■	■				■	■	
	9409047			■	■	■				■	■	■								■	■	■				■	■	
100-1000 µl long	9409048						■						■	■	■	■												
	9409049						■						■	■	■	■												

Compatibility of LLG-Pipette tips "economy" with Socorex Pipettes

Pipettor		Socorex											
Brand		Calibra digital 822 micro					Calibra digital 832 macro	Calibra digital 852 multi					
Model		variable/single channel					variable/single channel	variable/multi channel					
Volume Range		0.2-2 µl	1-10 µl	1-10 µl	2-20 µl	10-100 µl	20-200 µl	100-1000 µl	0.2-2 ml	1-10 ml	1-10 µl	10-100 µl	20-200 µl
Pipette tips "economy" line													
Volume	Cat-No.												
0.1-10 µl	9409099	■									■		
	9409045	■									■		
1-200 µl	9409046		■	■	■	■	■				■	■	
	6253577		■	■	■	■	■				■	■	
	9409047		■	■	■	■	■				■	■	
1-5 ml	6254364								■				
1-10 ml	9409098									■			

Compatibility of LLG-Pipette tips "economy" with HTL Pipettes

Pipettor		HTL																															
Brand	Labmate pro	Clinipet ⁺																															
Model	variable/single channel	fix/single channel																															
Volume Range	0.2-2 µl (red) 0.5-10 µl (red) 2-20 µl (yellow) 10-100 µl (yellow) 20-200 µl (yellow) 100-1000 µl (blue) 1000-5000 µl (white)	2 µl (yellow)	5 µl (yellow)	7 µl (yellow)	10 µl (yellow)	15 µl (yellow)	20 µl (yellow)	25 µl (yellow)	30 µl (yellow)	40 µl (yellow)	44.7 µl (yellow)	50 µl (yellow)	60 µl (yellow)	70 µl (yellow)	75 µl (yellow)	80 µl (yellow)	90 µl (yellow)	120 µl (yellow)	150 µl (yellow)	200 µl (yellow)	200 µl (blue)	220 µl (blue)	250 µl (blue)	300 µl (blue)	400 µl (blue)	450 µl (blue)	500 µl (blue)	600 µl (blue)	700 µl (blue)	750 µl (blue)	800 µl (blue)	900 µl (blue)	1000 µl (blue)
Pipette tips "economy" line																																	
Volume	Cat-No.																																
0.1-10 µl	9409099	■	■																														
	9409045	■	■																														
1-200 µl	9409046			■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
	6253577			■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
	9409047			■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
100-1000 µl long	9409048																				■	■	■	■	■	■	■	■	■	■	■	■	
	9409049																				■	■	■	■	■	■	■	■	■	■	■	■	■
1-5 ml	9409097																																■