

Compatibility of LLG-Tips Economy with LLG Pipettes

Model		variable/single channel							fix/single channel															
Volume Range / Cat. No.		0.1-2.5 µl / 9.280 000	0.5-10 µl / 9.280 001	2-20 µl / 9.280 002	10-100 µl / 9.280 003	20-200 µl / 9.280 004	100-1000 µl / 9.280 005	500-5000 µl / 9.280 006	2000-10000 µl / 9.280 007	5 µl / 9.277 515	10 µl / 9.277 516	20 µl / 9.277 517	25 µl / 9.277 518	50 µl / 9.277 519	100 µl / 9.277 520	200 µl / 9.277 521	250 µl / 9.277 522	500 µl / 9.277 523	1000 µl / 9.277 524	2000 µl / 9.277 525	2500 µl / 9.277 526	5000 µl / 9.277 527	10000 µl / 9.277 528	
Economy line Tips																								
Volume	Cat-No.																							
0.1-10 µl	9.409 099	■	■	■																				
	9.409 045	■	■	■																				
0.5-10/20 µl	6.252 446		■																					
	6.240 598		■																					
1-200 µl	9.409 046				■	■				■	■	■	■	■	■	■								
	6.253 577				■	■				■	■	■	■	■	■	■								
	9.409 047				■	■				■	■	■	■	■	■	■								
100-1000 µl long	9.409 048						■									■	■		■					
	9.409 049						■									■	■		■					
1-5 ml	9.409 097							■												■		■		
	6.254 364																				■	■		■
1-10 ml	9.409 098																						■	