

Technical Specification

PRODUCT : SILICA GEL 60 Å 35 -70 µm RS - For chromatography

CODE : P20000 **MOLECULAR WEIGHT :60,09**

CAS N° : 7631-86-9

FORMULA : SiO2

TEST	U.M.	SPECIFICATION
Particle size	µm	35 - 70
Particle size < 40µm	%	<= 15
Particle size > 63µm	%	<= 10
Loss on drying	%	<= 8.0
pH (5%w/w solution)	-	6.5 - 7.5

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in