

# Technical Specification

**PRODUCT : REAGENT TAN RS - For analysis**

**CODE : PS0327.**

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Colour	APHA	<= 10
Identification	-	Positive
Water (K.F.)	ppm	4500 - 5500
Acidity	ppm	<= 5

1L Content :  
- 495 mL Propanol-2  
- 500 mL Toluene  
- 5 mL HPLC water

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.  
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in