



	"A"	"B"	"C"	"D"	GROSS WEIGHT	NET WEIGHT
Safefast Classic 209	887	1045	427	1110	Kg. 180	Kg. 155
Safefast Classic 212	1192	1350	580	1470	Kg. 205	Kg. 175
Safefast Classic 218	1802	1960	885	2060	Kg. 280	Kg. 240

DESCRIZIONE  
**Safefast Classic**

CODICE DISEGNO: K6 00114 002\_COMM.    REVISIONE: 00    DEL: 21/07/10    DIS.: P.M.