



# Cloroformio

Synonyms:

Triclorometano

| CHCl <sub>3</sub> | Classification                  | Danger                                                                                                                                                                     |
|-------------------|---------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Molecular weight  | 119 ONU                         | 1888   |
| CAS               | 67-66-3 Transport hazard class: | 6.1                                                                                                                                                                        |
| EEC-N:            | 200-663-8 Packing               | III                                                                                                                                                                        |

Cloroformio

&gt; RE - Puro - Stabilizzato Etanolo

Grado Tecnico

|                          |                         |      |
|--------------------------|-------------------------|------|
| Description              | Clear colourless liquid | -    |
| Colour                   | <= 10                   | APHA |
| Identification (I.R.)    | Positive                | -    |
| Density at 20°C          | 1.478 - 1.488           | -    |
| Refractive index at 20°C | 1.4440 - 1.4460         | -    |
| Boiling point            | 61.0 - 61.5             | °C   |
| Acidity (HCl)            | <= 50                   | ppm  |



|                        |         |     |
|------------------------|---------|-----|
| Ethanol                | 0.6 - 1 | %   |
| Residue on evaporation | <= 20   | ppm |
| Water (K.F.)           | <= 0.03 | %   |
| Assay (GC)             | >= 99   | %   |

| ID                 | Size  | Packaging          |
|--------------------|-------|--------------------|
| 000000000000508320 | 1 l   | Bottiglia di vetro |
| 000000000000508321 | 5 l   | Tanica di plastica |
| 000000000000508322 | 200 l | Fusto di metallo   |
| P02402E48          | 25 l  | Fusto di metallo   |



CARLO ERBA Reagents S.r.l.  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



CARLO ERBA Reagents  
operates with a Certified  
Quality Management  
System

