



Sodium tartrate dihydraté

Synonyms:

L - (+) acide tartrique sel disodique ,
Tartrate disodique dihydraté

(CHOHCOONa) $\cdot$ 2.2H $_2$ O

Classification

Danger

Molecular weight **230.08** ONU

CAS **6106-24-7** Transport hazard
class:

EEC-N: **212-773-3** Packing

**Sodium tartrate
dihydraté**

> RPE - Pour analyse - ACS - Reag. Ph.Eur. -
Reag. USP

Grade Analytique

Description	White crystalline powder	-
Identification	Positive	-
pH sol. 5% at 25°C	7.0 - 9.0	-
Loss on drying (150°C)	15.61 - 15.71	%
Ammonium	<= 30	ppm
Chloride	<= 5	ppm



Sodium tartrate dihydraté > RPE - Pour analyse - ACS - Reag. Ph.Eur. - Reag. USP		Grade Analytique
Phosphate	<= 5	ppm
Insoluble in water	<= 50	ppm
Heavy metals (as Pb)	<= 5	ppm
Sulfate	<= 50	ppm
Ca	<= 100	ppm
Fe	<= 10	ppm
Assay (on dry)	99.0 - 101.0	%

ID	Size	Packaging
0000000000000483635	100 g	Flacon plastique
0000000000000483636	500 g	Flacon plastique
0000000000000483637	1 kg	Flacon plastique



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val de
 Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA Reagents
 operates with a Certified
 Quality Management
 System

