



Xilene, isomeri misti

Synonyms:

C8H10**Classification****Danger**

Molecular weight

106.17 ONU**1307**

CAS

Transport hazard
class:**3**

EEC-N:

905-588-0 Packing**III****Xilene, isomeri misti****> RE - Puro - Basso contenuto di benzene****Grado Tecnico**

Description	Clear liquid	-
Identification (I.R.)	Positive	-
Density at 20°C	0.862 - 0.872	-
Refractive index at 20°C	1.4917 - 1.5017	-
Residue on evaporation	<= 100	ppm
Water (K.F.)	<= 200	ppm



Benzene	<= 50	ppm
Total sulphur	<= 100	ppm
Assay (isomeric mixture)	>= 98.5	%

ID	Size	Packaging
000000000000392602	1 l	Bottiglia di vetro
000000000000392603	2.5 l	Bottiglia di vetro
000000000000392605	23 kg	Fusto di metallo
000000000000528251	5 l	Tanica di plastica
000000000000528252	200 l	Fusto di metallo
000000000000528254	10 l	Tanica di plastica