

**Bacinelle da laboratorio, acciaio 18/10****RSG**

Bacinelle da laboratorio, acciaio 18/10. Lucidate. Forma tonda con bordo e fondo piatto.



Code	Description	Packaging
LLG07605523	Bacinelle da laboratorio, acciaio 18/10, 14000 ml, Ø int. 400 mm, H 160 mm	1 pz
LLG09200518	Bacinella in acciaio inox, forma tonda, cap. 1000 ml Ø est. 165 mm h 75 mm	1 pz
LLG09200524	Bacinella in acciaio inox, forma tonda, cap. 2000 ml Ø est. 220 mm h 75 mm	1 pz
LLG09200526	Bacinella in acciaio inox, forma tonda, cap. 3000 ml Ø est. 240 mm h 95 mm	1 pz
LLG09200529	Bacinella in acciaio inox, forma tonda, cap. 4500 ml Ø est. 260 mm h 120 mm	1 pz
LLG09200533	Bacinella in acciaio inox, forma tonda, cap. 6000 ml Ø est. 305 mm h 120 mm	1 pz
LLG09200536	Bacinella in acciaio inox, forma tonda, cap. 8000 ml Ø est. 320 mm h 145 mm	1 pz
LLG09200541	Bacinelle da laboratorio, acciaio 18/10, Capacità 11000 ml, Ø int. 380 mm, Altezza 135 mm	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

