

**Anello salvagoccia, PP****Duran**

Anello salvagoccia, PP. Per bottiglie ISO. Autoclavabile fino a 140 °C (vapore).



Code	Description	Packaging
LLG09209078	Anello salvagoccia, PP, verde, Filettatura: GL 45	10 pz
LLG09209079	Anello salvagoccia, in PP, grigio, GL45	10 pz
LLG09209080	Anello salvagoccia, in PP, giallo, GL45	10 pz
LLG09209082	Anello salvagoccia, PP, GL 32, blu, resiste fino 140 °C, autoclavabile	10 pz
LLG09209084	Anello salvagoccia, PP, GL 45, blu, resiste fino 140 °C, autoclavabile	10 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

