

**Beuta Erlenmeyer, in vetro DURAN®****Duran**

Beuta Erlenmeyer, vetro DURAN®. Vetro DURAN®. Bordo diritto. Per uso con tappi Kapsenberg. Vetro tipo 1 / vetro neutro, in conformità con USP, EP e JP. Autoclavabile.





CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

