

**Cilindri graduati, PMP, Classe B, graduazione blu****Kartell**

Cilindri graduati, PMP, Classe B, graduazione blu. Forma alta, con becco, trasparente. In accordo con DIN 12681/ISO 6706. Calibrati a 20°C (TC, In). L'esposizione a temperature fino a 121°C (autoclavaggio) non causa il superamento permanente dei limiti di tolleranza. Si raccomanda di lavare a temperature al di sotto dei 60° per preservare graduazioni ed iscrizioni.



Code	Description	Packaging
LLG09274770	Cilindri graduati in PMP, 10 ml, grad. 2 ml, toll. $0.2 \pm \text{ml}$, h 140 mm, ATTENZIONE! ORDINE MINIMO 2 pz O MULTIPLI DI 2 pz	1 pz
LLG09274771	Cilindri graduati in PMP, 25 ml, grad. 5 ml, toll. $0.5 \pm \text{ml}$, h 195 mm	1 pz
LLG09274772	Cilindri graduati in PMP, 50 ml, grad. 10 ml, toll. $1.0 \pm \text{ml}$, h 200 mm	1 pz
LLG09274773	Cilindri graduati in PMP, 100 ml, grad. 10 ml, toll. $1.0 \pm \text{ml}$, h 250 mm	1 pz
LLG09274774	Cilindri graduati in PMP, 250 ml, grad. 20 ml, toll. $2.0 \pm \text{ml}$, h 315 mm	1 pz
LLG09274775	Cilindri graduati in PMP, 500 ml, grad. 50 ml, toll. $5.0 \pm \text{ml}$, h 360 mm	1 pz
LLG09274777	Cilindri graduati in PMP, 2000 ml, grad. 200 ml, toll. $20 \pm \text{ml}$, h 535 mm ATTENZIONE! ORDINE MINIMO 2 CONF. O MULTIPLI DI 2 CONF.	1 pz