




Provette, fondo tondo, vetro Soda-lime

Provette, fondo tondo, vetro Soda-lime. Vetro soda-lime. Spesso, robusto. Con bordo liscio e fondo tondo.



Code	Description	Packaging
LLG09400004	Provetta in vetro, 40x8 mm	200 pz
LLG09400007	Provetta in vetro, 70x10 mm	100 pz
LLG09400009	Provetta in vetro, 75x12 mm	100 pz
LLG09400011	Provetta in vetro, 85x15 mm	100 pz
LLG09400015	Provetta in vetro, 10x100 mm	100 pz
LLG09400017	Provetta in vetro, 100x12 mm	100 pz
LLG09400020	Provetta in vetro, 100x14 mm	100 pz
LLG09400022	Provetta in vetro, 100x16 mm	100 pz
LLG09400030	Provetta in vetro, 130x14 mm	100 pz
LLG09400040	Provetta in vetro, 160x16 mm	100 pz
LLG09400050	Provette in vetro, 180x18 mm	100 pz

