




Setacci per mulini a coltelli SM 100 / 200 / 300

Retsch

Setacci per mulini a coltelli SM 100 / 200 / 300.



Code	Description	Packaging
LLG06240026	Setacci per mulini SM 100/200/300, fori quadrati, luce netta 2,00 mm	1 pz
LLG09738815	Setacci per mulini SM 100/200/300, fori quadrati, luce netta 4,00 mm	1 pz
LLG09738816	Setacci per mulini SM 100/200/300, fori quadrati, luce netta 6,00 mm	1 pz
LLG09738817	Setacci per mulini SM 100/200/300, fori quadrati, luce netta 8,00 mm	1 pz
LLG09738818	Setacci per mulini SM 100/200/300, fori quadrati, luce netta 10,00 mm	1 pz
LLG09738854	Setacci per mulini SM 100/200/300, fori Conidur, luce netta 0,25 mm	1 pz
LLG09738855	Setacci per mulini SM 100/200/300, fori Conidur, luce netta 0,50 mm	1 pz
LLG09738856	Setacci per mulini SM 100/200/300, fori Conidur, luce netta 0,75 mm	1 pz
LLG09738857	Setacci per mulini SM 100/200/300, fori Conidur, luce netta 1,00 mm	1 pz
LLG09738858	Setacci per mulini SM 100/200/300, fori Conidur, luce netta 1,50 mm	1 pz
LLG09739055	Setacci per mulini SM 100/200/300, fori quadrati, luce netta 20,00 mm	1 pz

