

**Standard Indice di Rifrazione****Reagecon Diagnostics**

Standard certificati con valori espressi in Indice di Rifrazione (R.I.) unità nD.

- I risultati dei test sono accreditati ISO17025 con un'incertezza di misura fino a ± 0.00014 nD unità per tutti i valori di Indice di Rifrazione.
- I valori di Indice di Rifrazione Nominale sono presi dalle tabelle pubblicate ICUMSA
- Possono essere utilizzati con tutte le marche di rifrattometri
- Validità estesa di 1 anno
- Utile bottiglia contagocce di alta qualità



Code	Description	Packaging
LLG06258375	Standard indice rifrazione 1,381149 (± 0.00004) 15ml	1 pz
LLG06258376	Standard indice rifrazione 1,420087 (± 0.00004) 15ml	1 pz
LLG06265142	Standard indice rifrazione 1,40978 (± 0.00004), 15ml	1 pz
LLG06265877	Standard indice rifrazione 1,34325 (± 0.00004), 15ml	1 pz
LLG06265878	Standard indice rifrazione 1,37233 (± 0.00004), 15ml	1 pz
LLG06265879	Standard indice rifrazione 1,39032 (± 0.00004), 15ml	1 pz
LLG06265880	Standard indice rifrazione 1,44193 (± 0.00004), 15ml	1 pz
LLG07658877	Standard indice rifrazione 1,347824 (± 0.00004), 15ml	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

