



**Bottiglie per campioni di latte**

**Gerber Instruments**

Bottiglie per campioni di latte.



| Code        | Description                                                            | Packaging |
|-------------|------------------------------------------------------------------------|-----------|
| LLG06266500 | Bottiglia campionamento latte 50ml, in vetro, 4 marcature 10/20/3040mm | 1 pz      |
| LLG06266502 | Bottiglia campionamento latte 50ml, in plastica, tappo a scatto        | 1 pz      |



**CARLO ERBA Reagents S.r.l.**  
 VIA R. MERENDI, 22  
 20007 CORNAREDO (MI)  
 TEL. +39 02 93 99 190  
 FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

