

**Sonde per omogenizzatori ad Ultrasuoni Sonopuls, titanio****Bandelin**

Sonde per omogenizzatori ad Ultrasuoni Sonopuls, titanio. Sonde in titanio. Per collegamento a sonde pluristadio / sonde ausiliarie.



Code	Description	Packaging
LLG06251521	Sonda tipo MS 2.0, Ø 2,0mm, punta Micro, per 0,25-20ml, per Sonpuls mini20	1 pz
LLG06274618	Sonde per Sonopuls HD 4400, Tipo TS 413, Per SH 400 G, 100-750 ml, Ø 13mm	1 pz
LLG07647594	Sonde per SH 100 G SH 200 G, Per 20-200 20-900 ml, Per HD 4100/4200, Ø13mm	1 pz
LLG09650157	Sonda in titanio TS 102, Ø 2,5mm, per SH 100 G	1 pz
LLG09650158	Sonda in titanio TS 103, Ø 3mm, per SH 100 G	1 pz
LLG09650159	Sonda in titanio TS 104, Ø 4mm, per SH 100 G	1 pz
LLG09650161	Sonde x SH 100 G SH 200 G, x 5-75 10-100 10-350ml, x HD 4050/4100/4200 Ø6mm	1 pz
LLG09650165	Sonda tipo MS 1.5, Ø 1,5mm, punta Micro, per 0,1-10 ml, per Sonpuls mini20	1 pz
LLG09650166	Sonda tipo MS 2.5, Ø 2,5 mm, punta Micro, per 0,5-25 ml, per Sonpuls mini20	1 pz
LLG09650169	Sonde per Sonopuls HD 4400, Tipo TS 416, Per SH 400 G, 250-1000 ml, Ø 16mm	1 pz
LLG09650239	Sonda tipo VS 70 T, Ø 13 mm, puntale lungo, per HD 2070/2200 HD 3100/3200	1 pz
LLG09650240	Sonda tipo VS 200 T, Ø 25,0 mm, puntale lungo, per HD 2200/3200 HD 3400	1 pz
LLG09650244	Sonda tipo MS 72, Ø 2mm, Micropunta, per HD 2070/2200 HD 3100/3200	1 pz
LLG09650245	Sonda tipo MS 73, Ø 3mm, Micropunta, per HD 2070/2200 HD 3100/3200	1 pz
LLG09650246	Sonda tipo KE 76, Ø 6 mm, punta conica, per HD 2070/2200 HD 3100/3200	1 pz
LLG09650247	Sonda tipo TT 13 FZ, Ø 13mm, piastra titanio c/perforaz., per HD 2200/3200	1 pz
LLG09650250	Sonda tipo TT 13, Ø13mm, piastra in titanio, per HD 2070/2200 HD 3100/3200	1 pz





CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

