

**Anello salvagoccia per tappi a vite Premium, TpCh260****Duran**

Anello salvagoccia per tappi a vite Premium, TpCh260. Anello salvagoccia incolore TpCh260, altezza 4 mm. Resistente a temperature da -196°C a +260°C.





CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

