



Conduttimetri P 2/30

Evoqua Water Technologies

Conduttimetri P 2/30.  
Con display analogico e attacco filettato. Con connessione a vite 3/4" e cella integrata (c=0.2).



Code	Description	Packaging
LLG09914516	Conduktivimetri P 2/30, 0,5 ... 20 $\mu\text{S}/\text{cm}$ , temp. 0 ... 50 °C, dim. (LxPxH) 106 x 62 x 126 mm	1 pz



**CARLO ERBA Reagents S.r.l.**  
 VIA R. MERENDI, 22  
 20007 CORNAREDO (MI)  
 TEL. +39 02 93 99 190  
 FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

