



Accessori per Data Logger PDF EBI 300 / EBI 310

EBRO

Accessori per Data Logger PDF EBI 300 / EBI 310



Code	Description	Packaging
LLG06263842	Sonda di ricambio per EBI 310 TH, tipo TPH 500	1 pz
LLG06281569	Filtro-Sinterizzato-Acciaio inox	1 pz
LLG09908099	Sonda di temperatura tipo TPX 120	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

