



**Accessori per Viscosimetri B-ONE TOUCH / FIRST TOUCH**

**Lamy Rheology**

Accessori per viscosimetri B-ONE TOUCH / FIRST TOUCH.



| Code        | Description                                                                 | Packaging |
|-------------|-----------------------------------------------------------------------------|-----------|
| LLG04659148 | Accessori per Viscosimetri B-ONE TOUCH / FIRST TOUCH, Tipo Axis R 1-6       | 1 pz      |
| LLG04659149 | Accessori per Viscosimetri, Tipo R-1, Disco, Ø 56,26mm, 50-600 000 mPas     | 1 pz      |
| LLG04659150 | Accessori per Viscosimetri, Tipo R-2, Disco, Ø 46,93mm, 100-2 400 000 mPas  | 1 pz      |
| LLG04659151 | Accessori per Viscosimetri, Tipo R-3, Disco, Ø 34,69mm, 100-6 000 000 mPas  | 1 pz      |
| LLG04659152 | Accessori per Viscosimetri, Tipo R-4, Disco, Ø 27,30mm, 200-12 000 000 mPas | 1 pz      |
| LLG04659153 | Accessori per Viscosimetri, Tipo R-5, Disco, Ø 21,14mm, 400-24 000 000 mPas | 1 pz      |
| LLG04659154 | Accessori per Viscosimetri, Tipo R-6, Disco, Ø14,62mm, 1000-60 000 000 mPas | 1 pz      |
| LLG04659155 | Accessori per Viscosimetri, Tipo Axis R-7, Ø 3,20mm, 4000-240 000 000 mPas  | 1 pz      |
| LLG04659156 | Accessori per Viscosimetri, Tipo Axis L-R, Asse adattamento                 | 1 pz      |
| LLG04659157 | Accessori per Viscosimetri, Tipo L-1, asta, Ø 18,80 mm, Lungh. 65,1 mm, 15- | 1 pz      |
| LLG04659158 | Accessori per Viscosimetri, Tipo L-2, asta, Ø 18,72 mm, Lungh. 6,86 mm, 50- | 1 pz      |
| LLG04659159 | Accessori per Viscosimetri, Tipo L-3, asta, Ø 12,60 mm, Lungh. 1,78 mm, 200 | 1 pz      |
| LLG04659160 | Accessori per Viscosimetri, Tipo L-4, asta, Ø 3,20 mm, Lungh. 31 mm, 1000-2 | 1 pz      |
| LLG04659161 | Accessori per Viscosimetri, Tipo L 1-4, asta                                | 1 pz      |





**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

